

CP♥

CARDIO PULMONARY  
VENTILATION



What if you could  
enhance your  
CPR performance?

## DEALING WITH CARDIO-RESPIRATORY ARREST:

### Adequate ventilation is needed to preserve circulation

Ventilation can adversely affect circulation while delivering chest compressions (CCs). Its role is, however, crucial when cardiopulmonary resuscitation (CPR) is prolonged. To meet CPR requirements, in line with international guidelines, ventilation should include high-quality CC and optimal oxygenation.<sup>1,2</sup>

The need for ventilation is low when CPR is prolonged, whereas optimal oxygenation is critical. The use of tracheal intubation, well-established as the "gold standard" of care in emergency, is currently being debated. Ventilation with mask interfaces is recommended, as long as the risks associated with bag valve ventilation are avoided.<sup>3,4</sup>

## AIR LIQUIDE MEDICAL SYSTEMS' COMMITMENT

For 40 years, Air Liquide Medical Systems has committed its expertise to healthcare professionals and patients to create and provide innovative medical devices, notably in respiratory care (ventilators and patient interfaces for intensive care, emergency, transport and patients' home).


To this end, Air Liquide Medical Systems developed the CPV (Cardio Pulmonary Ventilation) solution. The solution is dedicated to managing cardiac arrest and can be used with the well-established emergency transport ventilator Monnal T60. CPV offers technical innovations for each stage of cardiopulmonary resuscitation (CPR) that are specifically designed for rescuers. The solution can be used by mobile intensive care units, during inter- and intra-hospital transports, and by specialized services.

1. European Resuscitation Council Guidelines for Resuscitation, 2010 ; 2. American Heart Association (AHA) Guidelines for CPR and ECC, 2010 ; 3. Hasegawa K et al. JAMA 2013; 309 : 257–266 ; 4. Segal N et al. Resuscitation 2015; 86 : 62–66.



# CPV, THE SOLUTION FOR A SYNCHRONIZED VENTILATION

Developed with professional rescuers to meet their expectations



ONE TOUCH to start  
LIMITS the risks associated with bag valve ventilation  
ENHANCES circulation  
LIMITS alarms  
GUIDES users to CPR

## What is it?

- ♥ The CPV (Cardio Pulmonary Ventilation) mode is an innovative solution in cardiac arrest management. It can be provided with the emergency transport ventilator Monnal T60
- ♥ Its goals are to easily perform ventilation and improve circulation at the same time
- ♥ It facilitates CPR management

## How to use it & how it works?

- ♥ The single button, displayed on the home screen, immediately launches the CPV mode under preset settings in line with international guidelines

These settings can be modified according to clinical requirements, when necessary

- ♥ The CPV solution assists users throughout the entire process of CPR, from the early stages of CPR, such as ventilation with mask interfaces (to avoid the risks associated with bag valve ventilation), until the more demanding technical stages, such as tracheal intubation
- ♥ The CPV solution ensures that the CC and return of spontaneous circulation (ROSC) phases switch when necessary, and so facilitates ventilation management

**CP** ♥  
CARDIO PULMONARY  
VENTILATION

# CPV, PROVIDING SYNCHRONIZED VENTILATION

## Just a single touch to combine ventilation with h

### Real-time performance feedback on CC quality to guide you in the management of CPR

---

- ♥ CC quality is known to affect patient survival. The use of a real-time monitoring system to report on indicators of external cardiac massage during CPR is an essential step forward
- ♥ The CPV solution, which is available on the emergency transport ventilator Monnal T60, offers immediate assistance to deliver high quality CPR
- ♥ Four essential parameters are detected in real time:
  - fCC, chest compression frequency
  - no CC, cumulative time during which no CCs were performed
  - P-P, indicator of CC strength
  - etCO<sub>2</sub>, concentration of exhaled CO<sub>2</sub>

#### CPV, the solution for a synchronized ventilation:

- detects CCs, monitors essential CC parameters to obtain immediate feedback on the quality of CPR performance
- synchronizes ventilation with CCs for optimal circulation

# VENTILATION FOR CPR

## with high-quality chest compressions

### Ventilation is actively synchronized with CCs for optimal circulation

- ♥ The CPV mode includes a specific synchronization algorithm for high and low pressures, which, for the first time, magnifies intrathoracic positive pressure (ejection) during CCs and negative intrathoracic pressure (venous return) generated by chest recoil
- ♥ The CPV solution guarantees sufficient oxygenation without interrupting external cardiac massage (ECM) and maintains sufficient alveolar ventilation

- CPV helps you focus your attention and efforts on patient outcomes
- CPV aids you in performing CPR within the recommended range of international guidelines:
  - You are given a continuous feedback loop to evaluate the quality of your actions and are guided throughout the entire process to improve CPR in real time
  - The CPV synchronized ventilation is designed to offer better oxygenation and improve patient hemodynamics



# CPV, PROVIDING SYNCHRONIZED VENTILATION

## Essential CC indicators to help you when performing

### Frequency (fCC)

- Critical parameter reflecting CPR efficacy
- Helps target and maintain optimal CC frequency

#### YOUR SUPPORT

On screen: frequency of chest compressions per minute  
Goal: 100/min



### Exhaled CO<sub>2</sub> monitoring (etCO<sub>2</sub>)

- Indicator of circulation function

#### YOUR SUPPORT

Detects correct intubation  
and CPR efficacy  
ROSC indicator

# RELATION FOR CPR: IN PRACTICE

## performing CPR

### Detection of compression interruptions (no CC)

- Indicator of cumulative time during which no CCs were performed

#### YOUR SUPPORT

Encourages the rescuer not to interrupt CCs

### Strength (P-P)

- Indicator of CC efficacy (indirectly shows the pressure transmitted to the thorax)
- Ensures regularity of compression amplitude

#### YOUR SUPPORT

Encourages the rescuer to maintain the pressure intensity of compressions or ask a colleague to take over

### Insufflated tidal volume (VTi)

- Monitoring of insufflated tidal volume while performing continuous chest compressions

### Ventilation settings zone

- Two levels of pressure-controlled ventilation synchronized with CCs
- Initial preset parameters in line with international guidelines:
  - FiO<sub>2</sub>: 100%
  - Frequency: 10 cycles/min
  - Duration of insufflation: 1 seconde

## Contact

Air Liquide Medical Systems  
Parc de Haute Technologie  
6, rue Georges Besse  
92182 Antony Cedex, France  
Tel. : +33 (0)1 40 96 66 00  
Fax : +33 (0)1 40 96 67 00  
[infoCPV@airliquide.com](mailto:infoCPV@airliquide.com)



Air Liquide Healthcare is a world leader in medical gases, home healthcare, hygiene products and healthcare specialty ingredients. It aims to provide customers in the continuum of care from hospital to home with medical products, specialty ingredients and services that contribute to protecting vulnerable lives.