

AInest mC1



Light and Intuitive


 **Air Liquide**
HEALTHCARE

Alnest MC1

Feel better everyday

Minimal contact nasal cushion

Comfortable soft silicone cushion
that fits under the nose
Easy to assemble thanks to a visual
and tactile mistake proof system
2 sizes (S and M/L)




Minimalist headgear with Quick Adjust function

Designed to help the patient
to adjust the mask thanks
to the patented landmarks system
Soft, light & seamless fabric

Easy to handle and set

Easy-peel function for an
easier setting, even for
patient with low dexterity




Silent airflow

A minimal leak noise
for a peaceful sleep



Flexible tube

For a complete freedom
of movement





Enhanced patient experience

With its 42 grams and compact design, MC1 fits for patients needing **discretion, comfort and simplicity.**

MC1 is one of the **lightest** nasal minimum contact masks on the market (only 42 g).¹



Discreet

The leak has been designed to **minimize the noise** and the diffusion of exhaled air, promoting a **better sleep quality** for both patient and partner.



Comfort

Thanks to its soft sub-nostril cushion, the MC1 mask sits comfortably **under the nose**, reducing red marks on the skin and completely removing the pressure on the nasal bridge. Its **open-view design** offers the patient a clear field of vision, allowing him to **wear glasses, read and watch TV**. This can significantly improve the quality of life for patients requiring home ventilation.



Easy to use

The landmarks make the mask fitting effortless for both patient and homecare provider. **For the patient**, a correct headgear fit is essential both to **maintain the mask seal** and to **reduce pressure points**. The landmarks thus promote **patient adherence** to treatment as a well-fitted mask is more likely to be worn regularly. **For the provider**, this innovation **saves time** during the initial installation, **reduces the number of follow-up visits** and enables the **remote patient assistance**.

Alnest MC1



Watch
the video



Mask references

Product reference	Description	Size	Distributor code
472113	Alnest MC1	S	KM229157
472114	Alnest MC1	M-L	KM229158


Manufactured by

AIR LIQUIDE MEDICAL SYSTEMS S.r.l
Via Bisceglie, 66 - 20152 Milano (MI) - ITALY
PLANT: Via dei Prati, 62 - 25073 Bovezzo (BS) - ITALY
Tel. +39 030.201.5911
Fax +39 030.209.83.29

Distributor contact

AIR LIQUIDE MEDICAL SYSTEMS S.A.
Parc de Haute Technologie
6, rue Georges Besse
92182 Antony Cedex
France

medicaldevice.airliquide.com

 **Hotline ALMS: +33 179 5170 01**

 **Air Liquide**
HEALTHCARE

ALNEST MC1 is a reusable nasal mask with calibrated exhalation orifices (Vented), intended to be used at home or in a hospital setting by a single adult patient (weighing more than 30 kg) who has been prescribed for a non-invasive positive pressure ventilation therapy (NPPV), for example a continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

Medical devices class IIa CE0051. Please read the user manual.

¹Comparative Analysis of the Weight of Leading Minimal Contact Masks for CPAP and BiPAP Delivery on the Market_Marius_L_2023