


Alnest

N1 Silent
N1



Silent and intuitive

UNIVEO Δ headgear

The first innovative headgear designed to help the patient

Quick Adjust, the Univeo Delta headgear scaling feature (patent filed)



Landmarks on the lower straps

Quick Adjust feature
with landmarks
Easy to handle and set
Soft & seamless fabric
Two-tone

- ▶ Facilitates the installation and adjustment of the mask
- ▶ Promotes the patient compliance to treatment

A generation of innovative masks

Alnest N1 Silent N1

Univeo Delta headgear

With **Quick Adjust** scaling function

A headgear designed to help patients

Available in S, M, L



Frame

Minimalist and light
Optimal coupling ergonomics

Quiet Flow

The function that ensures optimal silence
Only 12.3 ± 3.2 dB(A) sound pressure level

N1 Silent

Secure Clip

The mistake-proof system

Cushion

Optimal sealing and comfort
Available in S, M and L sizes

N1

Also available with standard leak (without Quiet Flow)

Ergonomic hooks

Quick assembly and disassembly of the headgear



Made in Italy

Alnest N1 Silent N1

3 headgear
sizes

3 cushion sizes

1 frame

Product Reference	Description	Size	Distributor code
472367	ALNEST N1 Silent	S	KM229116
472368	ALNEST N1 Silent	M	KM229117
472369	ALNEST N1 Silent	L	KM229118
472486	ALNEST N1	S	KM229148
472487	ALNEST N1	M	KM229149
472488	ALNEST N1	L	KM229150

Spare parts:

472478	N1 & N1 Silent cushion	S	KM229122
472479	N1 & N1 Silent cushion	M	KM229123
472480	N1 & N1 Silent cushion	L	KM229124
472481	Univeo Delta headgear	S	KM229125
472482	Univeo Delta headgear	M	KM229126
472485	Univeo Delta headgear	L	KM229131
472483	Swivel connector	-	KM229130


Manufactured by

AIR LIQUIDE MEDICAL SYSTEMS S.r.l
Via Bisceglie, 66 - 20152 Milano (MI) - ITALY
PLANT: Via dei Prati, 62 - 25073 Bovezzo (BS) - ITALY
Tel: +39 030.201.5911
Fax +39 030.209.83.29

Distributor contact

AIR LIQUIDE MEDICAL SYSTEMS S.A.
Parc de Haute Technologie
6, rue Georges Besse
92182 Antony Cedex
France.

medicaldevice.airliquide.com
Mail: dl-frtaema-src-export@airliquide.com

 **Hotline ALMS: +33 179 51 70 01**



 **Watch
the video**

 **Air Liquide**
HEALTHCARE

ALNEST N1 and ALNEST N1 Silent are reusable nasal masks with calibrated exhalation orifices (Vented), intended to be used at home or in a hospital setting by multiple adult patients (weighing more than 30 kg) who have been prescribed for a non-invasive positive pressure ventilation therapy (NPPV), for example a continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

*ALNEST N1 Silent and ALNEST N1 are class IIa medical devices CE0051. Please read the respective user manual.