





UNIVEO™ Headgear

the silicone bumper

• 5 connection points for a peripheral stability and even distribution of pressure points

Self-regulation of the inclination

of the frontal branch thanks to



Thin and flexible frame

- Provides an unobstructed field of vision
- Optimal contact of the cushion on the patient's face thanks to the flexibility of the frontal branch

Quiet leak

• Provides silent CO2 washout for a comfortable sleep experience for both the patient and his/her partner



360° rotating elbow

 Allows for an enhanced freedom of movement to minimize any disruptions during the patient's sleep

RespireoTM SOFT Facial is also available in Non Vented version

Scan here for

more information



Air Liquide

HEALTHCARE

Cushion SOFT

- Developed in close partnership with experts
- Ensures a balance between flexibility and firmness while providing reliable sealing at varying pressures and comfort for the patient

available S, M and L



• For a quick assembly and disassembly of the UNIVEOTM headgear

Product Reference	Description	Size	Distributor code
472548	Respireo™ SOFT Facial	S	KM229103
472549	Respireo™ SOFT Facial	M	KM229104
472550	Respireo™ SOFT Facial	L	KM229105
472551	Respireo™ SOFT Facial Non Vented	S	KM229119
472552	Respireo™ SOFT Facial Non Vented	M	KM229120
472553	Respireo™ SOFT Facial Non Vented	L	KM229121
472557	Respireo™ SOFT Facial Annual pack	S	KM229113
472558	Respireo™ SOFT Facial Annual pack	M	KM229114
472559	Respireo™ SOFT Facial Annual pack	L	KM229115

Spare parts:

472471	UNIVEO™ Headgear	-	KM227700		
472633	Facial cushion	S	KM229110		
472634	Facial cushion	М	KM229111		
472635	Facial cushion	L	KM229112		
472483	Grey swivel fitting*	-	KM229130		

Manufactured by

AIR LIQUIDE MEDICAL SYSTEMS S.r.I Via Bisceglie, 66 - 20152 Milano (MI) - ITALY PLANT: Via dei Prati, 62 - 25073 Bovezzo (BS) - ITALY Tel +39 030 201 5911 Fax +39 030.209.83.29

Distributor Contact

AIR LIQUIDE MEDICAL SYSTEMS S.A. Parc de Haute Technologie 6, rue Georges Besse 92182 Antony Cedex France

medicaldevice.airliquide.com

Respireo™SOFT Facial is a reusable oronasal mask with calibrated exhalation orifices (Vented), intended to be used at home or in a hospital setting by multiple adult patients (weighing more than 30 kg) who have been prescribed for a non-invasive positive pressure ventilation therapy (NPPV), for example a continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

Respireo™ SOFT Facial Non Vented is a reusable oronasal mask without calibrated exhalation holes intended to be used at home or in a hospital setting by multiple adult patients (weighing more than 30 kg) under Non invasive Positive-Lorent None (weighing more than 20 kg) under None (weighing more than 30 kg) under None (weighing moPressure Ventilation (NPPV).

* This spare part is for the Vented version.

Medical devices class IIa CE0051. Please read the user manual.

