

Respireo

Primo P



Simple & Ergonomic pillow mask



We are all different. The Respireo Primo series is a range of masks for ventilation that has been **specifically designed to fit each need of the patients and adapt to each morphology.**

Simple & Ergonomic, **Respireo Primo P** is a pillow mask combining silence and ease of use in order to quickly win the compliance of the patients.

Available in vented and non vented versions, it is the reference mask for all patients willing to keep their liberty and autonomy in their daily activities.

It is also recommended for patients inconvenienced by forehead support, contact on the nasal bridge or suffering from claustrophobia.



Respireo

Primo P



Adaptable tubing

- Ball joint system for a rotation at 360°
- Adjustable fixation: over the head, on the lateral sides or free downwards
- Clear visual field



Ergonomic cushion

- Airtight and fitting the nostrils thanks to the large sides and the "bellows" technology
- Easy assembly of the cushion with a visual and tactile poka-yoke

Stable headgear

- Designed to ensure settings, comfort and stability without any risk of sliding back of the head
- Two levels of adjustments



Silent and diffuse leak



Respireo Primo P P is also available in **non vented version** for patients ventilated with valve ventilators.



Codes	Désignation	Size
KM219900	Respireo Primo P (manufacturer reference 472105*)	S, M, L
KM220000	Respireo Primo P Non Vented (manufacturer reference 472106*)	S, M, L
KM220700	Cushion Respireo (manufacturer reference 472210)	S
KM220800	Cushion Respireo Primo P (manufacturer reference 472211)	M
KM220900	Cushion Respireo Primo P (manufacturer reference 472212)	L



Please read user manual.

(*) Medical device class IIa CE0051 manufactured by Air Liquide Medical Systems S.r.l.
Distributed by Air Liquide Medical Systems S.A.